

## PB 735

(PB735L-4C-90200-S-DBFK-C4, article no.: 83751)

TECHNICAL SPECIFICATION		
<i>Load</i> [kg]	Max. Patient weight	180
	Safe Working Load	220
	Lifting pole load	80
<i>Dimensions</i> [cm]	External [W x L]	approx. 107 x 206
	Mattress support platform [W x L]	90 x 200
	Under bed clearance	approx. 10,5
<i>Adjustment</i>	Height [cm]	approx. 31,0 ÷ 81,0
	Angle – back rest section [°]	approx. 0 ÷ 70
	Angle – upper leg rest section [°]	approx. 0 ÷ 32
	Angle – lower leg rest section [°]	approx. 0 ÷ 20
<i>Mattress</i>	Angle – anti-Trendelenburg position [°]	approx. 0 ÷ 12
	Dimensions [HxWxL]	20x90x200
<i>Drive system</i>	Volume density of foamed material	35 - 50 kg/m <sup>3</sup>
	Drive lifting unit	Linak LA40
	Mattress support platform drive unit	Linak LA27
	Control box	Linak C061
	Hand switch	HL754
	Input voltage	230V AC, 50Hz
	Operating time (duty cycle)	2 min on / 18 min off
	Ingress protection	IP X4
	Electrical appliance class	II
	Noise level	< 65 dB (A)
<i>Weight</i> [kg]	Total	143,0
	Mattress support platform (head section)	50,2
	Mattress support platform (foot section)	10,5
	End board (1)	41,0
	Scissor lift with drive unit	1,9
	Mattress support platform drive unit (1)	1,8 / 1,5
	Hand switch / control box	0,3 / 0,9
	Side rails (set)	17,2
	Lifting pole	5,3
<i>Conditions operating</i>	Temperature [°C]	+10 ÷ +40
	Humidity [%]	30 ÷ 75
<i>Conditions storage</i>	Temperature [°C]	+5 ÷ +45
	Humidity [%]	30 ÷ 75



<i>Programmable functions</i>	The height adjustment intermediate stop	approx. 38 cm
	The back rest section adjustment intermediate stop	approx. 30°

According to the condition of **5<sup>th</sup> March 2021**